User Guide/Manuel d'utilisation/Manuale d'uso/Benutzerhandbuch/
Guía del usuario/ユーザー マニュアル / 用户指南 / 사용자 설명서

Millex®-GP (33 mm)

Sterilizing Filter Unit

Unité de filtration stérilisante
Unità di filtrazione sterilizzante
Sterilfiltereinheit
Unidad de filtro esterilizante
滅菌フィルター ユニット

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SLGP 033 RS (50/Box/Boîte/Conf./Packg./por caja/箱 /箱 / 상자)

SLGP 033 RB (250/Box/Boîte/Conf./Packg./por caja/箱 /箱 / 상자)

MILLIPORE

- Single use only
- Sterile
- Non-pyrogenic
- Not for use as an in-line filter for intravenous administration.

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Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore Corporation. Millipore Corporation assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore Corporation be liable for incidental or consequential damage in connection with or arising from the use of this manual.

Copyright 2005, Millipore Corporation. All rights reserved. Printed in Ireland.

Millipore, Millipore Express, and Millex are registered trademarks of Millipore Corporation.

Freon is a trademark of E. I. du Pont de Nemours and Company, Inc.

Luer-Lok is a trademark of Becton, Dickinson and Company.

Alconox is a registered trademark of Alconox, Inc.

Lubrol is a registered trademark of Imperial Chemical Industries, Ltd.

Nonidet is a trademark of Royal Dutch/Shell Group.

Tween is a registered trademark of Atlas Powder Company.

P36359, Rev C, 11/05

This product is manufactured by:
Millipore Ireland BV,
 Carrigtwohill,
 Co. Cork

Language	Page No.
English	1
Français	10
Italiano	19
Deutsch	28
Español	37
日本語	46
中文	51
한국어	56

Introduction

This document provides compatibility information, operating steps, and specifications for the Millipore Express® PES family of sterile Millex® filter units. Millex filter unit's bi-directional support of the filter membrane enables users to filter aqueous solutions in either direction, forward (from the syringe into the container) or backward (from the container into the syringe). The Millex-GP filter unit removes microorganisms, particles, precipitates, and undissolved powders larger than 0.22 micron (μm). Typical research laboratory applications include the sterile filtration of protein solutions, tissue culture media, additives, buffers, and water. These single-use filter units consist of a membrane filter sealed in an acrylic housing. They are non-pyrogenic and non-toxic.

PLEASE NOTE:

Outside of the U.S. and Japan, Millex-GP filter units (Cat Nos. SLGP 033 RS and SLGP 033 RB) may also be used for direct patient care and pharmacy admixture applications. In the hospital pharmacy, Millex-GP filter units can be used to sterile filter small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization and particulate removal from epidural and other liquid anesthetics as well as from irrigation solutions used in ophthalmic, otic, and other surgical procedures.

Within the U.S. and Japan, the approved Millex-GP filter unit for direct patient care and pharmacy admixture applications is Cat. No. SLGP M33 RS.

Chemical Compatibility

The Millex-GP filter unit is compatible with most aqueous solutions. Based on information from technical publications, materials suppliers, laboratory tests, and field evaluations, Millipore believes that the agents listed in the following chart are safe to use with Millex filter units. However, because of the effects of variability in temperature, concentrations, duration of exposure, and other factors outside of our control, Millipore does not provide or imply a warranty with respect to this information. Agents that are not listed should be tested prior to use.

Chemical

Acetic Acid (aqueous solution)	Hydrochloric Acid (1N solution)	Pentane
Alconox® detergent 1%	Hydrochloric Acid (6N solution)	Petroleum Ether
Ammonium Hydroxide	Hydrochloric Acid (concentrated)	Phenol (aqueous solution)
Ammonium Sulfate (saturated)	Hydrofluoric Acid	Potassium Hydroxide (3N solution)
Boric Acid (aqueous solution)	Hydrogen	Silicone Oils
CHAPS (aqueous solution)	Hydrogen Peroxide (3% solution)	Sodium Carbonate (aqueous solution)
Diethyl Pyrocarbonate (0.2% solution)	Hydrogen Peroxide (30% solution)	Sodium Chloride (2M solution)
Ethylene Glycol	Hydrogen Peroxide (90% solution)	Sodium Dodecyl Sulfate
Formaldehyde	HYPO (dilute solution)	Sodium Hydroxide (3N solution)
Formic Acid (50% solution)	Kerosene	Sodium Hydroxide (concentrated)
Freon® solvent (TF or PCA)	Lactic Acid (50% solution)	Sulfuric Acid (6N solution)
Gasoline	Lubrol® PX emulsifier	Trichloroacetic Acid (aqueous solution)
Glycerine (glycerol)	Mercaptoethanol (0.1M solution)	Tween® 20 emulsifier
Guanidine HCl (6M solution)	Mineral Spirits	Urea (8M solution)
Guanidine Thiocyanate (5M solution)	Nitric Acid (6N solution)	Water (Brine - Sea Water)
Guanidinium HCL (6M solution)	Nitric Acid (concentrated)	Water (de-ionized)
Helium	Nitrogen	
Hexane	Nonidet™ -P 40 surfactant	
	Paraldehyde	

Active Drug Compounds

Drug	Ref. No.	Merck® Index 11th Edition Page No.	Express (PES) Membrane Compatibility	Comments
Aminophylline	477	76	OK	water soluble
Ampicillin	621	93	OK	water soluble
Aspartame	861	132	OK	water soluble
Bleomycins	1324	201	OK	water soluble
Caffeine	1635	248	OK	water soluble
Cefazolin	1925	294	OK	water soluble
Cefoxitin	1938	297	OK	water soluble
Cephalothin	1978	305	OK	water soluble
Cisplatin	2319	361	OK	water soluble
Colistin	2475	387	OK	suspension + surfactant
Cytarabine	2790	437	OK	water soluble
Dactinomycin	2804	441	OK	water soluble
Daunorubicin	2825	445	OK	water soluble
Dexamethasone	2922	463	OK	5% alcohol
Diazepam	2977	472	OK	40% alcohol
Digoxin	3150	499	OK	50% alcohol
Dobutamine	3396	535	OK	water soluble
Dopamine	3415	538	OK	water soluble
Doxorubicin	3428	540	OK	water soluble
Ergonovine	3600	573	OK	water soluble
Etoposide	3842	610	OK	30% alcohol
Factor III	3873	616	OK	water soluble
Factor IX	3874	616	OK	water soluble
Fentanyl	3944	628	OK	water soluble
Fluorouracil	4109	654	OK	water soluble

Active Drug Compounds, continued

Drug	Ref. No.	Merck Index 11th Edition Page No.	Express (PES) Membrane Compatibility	Comments
Folic Acid	4140	660	OK	water soluble
Furosemide	4221	674	OK	water soluble
Gentamicin	4284	686	OK	water soluble
Hemin	4563	733	OK	water soluble
Heparin	4571	735	OK	water soluble
Hydrocortisone 21-glycol			OK	100% propylene glycol
Sodium Succinate	4713	758	OK	water soluble
Immunoglobulins	4837	780	OK	water soluble
Insulin	4887	789	OK	water soluble
Isoproterenol	5105	821	OK	water soluble
Lidocaine	5359	863	OK	water soluble
Mannitol	5629	901	OK	water soluble
Metronidazole	6079	968	OK	water soluble
Mitoguazone	6131	978	OK	water soluble
Mitomycins	6133	979	OK	water soluble
Mitoxantrone	6135	979	OK	water soluble
Moxalactam	6201	991	OK	water soluble
Nitroglycerin	6528	1045	OK	water soluble
Norepinephrine	6612	1058	OK	water soluble
Penicillin G potassium	7041	1123	OK	water soluble
Phenobarbital	7201	1149	OK	water soluble
Piperacillin	7430	1184	OK	water soluble
Plicamycin	7510	1198	OK	water soluble
Prochlorperazine	7768	1231	OK	water soluble
Protamines	7898	1253	OK	water soluble

Active Drug Compounds, *continued*

Drug	Ref. No.	Merck Index 11th Edition Page No.	Express (PES) Membrane Compatibility	Comments
Streptokinase	8784	1390	OK	water soluble
Sulfamethomidine	8888	1407	OK	50% alcohol
Tobramycin	9413	1494	OK	water soluble
Trimethoprim	9624	1528	OK	water soluble
Urokinase	9799	1555	OK	water soluble
Vidarabine	9881	1569	OK	water soluble
Vinblastine	9887	1570	OK	water soluble
Vincristine	9891	1571	OK	water soluble

How to Use the Millex Sterilizing Filter Unit

This section lists warnings and cautions and provides steps to use the Millex sterilizing filter unit.

WARNINGS:

- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous fluid administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 cc because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.

CAUTIONS:

- Do not use the Millex filter unit to filter fluids at temperatures above 45 °C (113 °F).
- Do not use the Millex filter unit to filter emulsions or suspensions because it was not designed for that purpose.
- Do not use it to filter solutions containing 5 milligrams (mg) or less of active drug materials unless binding studies have been performed.
- Do not use the same Millex filter unit to filter solutions in both directions.
- Do not re-sterilize or reuse the Millex filter unit. It is a single-use only device.

Procedure for Using the Millex Sterilizing Filter Unit

To view the steps for using a Millex sterilizing filter unit, see the “Overview of Procedure” section at the back of this user guide.

Specifications for 0.22 μm Millex-GP Filter Unit

Materials	Hydrophilic polyethersulfone membrane; acrylic housing	
Dimensions	inlet to outlet	26 mm (1.02 in.)
	diameter	33 mm (1.30 in.)
	filtration area	4.52 cm ² (0.70 in ²)
	pore size	0.22 μm
Temperature limit	45 °C (113 °F) maximum	
Housing burst at 25 °C	150 psig (10 bar) inlet maximum	
Filtration volume	10 mL to 200 mL	
Hold-up volume	≤ 0.1 mL after air purge	
Connections	Female Luer-Lok® inlet; male Luer-slip outlet	
Sterilization method	Gamma irradiation	
Flow rate at 25 °C	Typical initial flow rate for water at 30 psig: 175 mL/min	

Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., see your Millipore catalogue for the phone number of the office nearest you or go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

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